


- 1 Upon receipt of your battery, it is recommended to bring the battery to complete charge for 12 hours prior to first use.




- 2 To charge the battery, plug the AC adaptor that comes with your Cholestech LDX™ analyzer into a wall outlet.



- 3 Connect the battery to the AC adaptor using the appropriate port on the AC adaptor and battery.



- 4 Flashing lights will indicate the battery is charging.



- 5 When the lights are completely full (5 LEDs on) the battery is fully charged.



When the charge is 100% the LEDs will turn off. Remove from the charger when all LEDs are off.

- 6 The battery is ready for use.



LED quick reference guide

In order to check the status of the battery charge level, press the  button. If the battery is being charged, the LEDs will display automatically.

- 1 1 red LED on and flashing indicates the battery must be charged – a test on the LDX might not complete



- 2 1 LED on not flashing indicates the battery level is <20% and should be charged soon



- 3 2 LEDs on indicates the battery level is <40%



- 4 3 LEDs on indicates the battery level is <60%



- 5 4 LEDs on indicates the battery level is <80%



- 6 5 LEDs on indicates the battery level is between 80% and 100%



- 7 LEDs flashing, while plugged in, indicates the battery is charging and the approximate charge level it is currently at (1-5 LEDs)



When the charge is 100% the LEDs will turn off. Remove from the charger when all LEDs are off.

Battery specifications

- The battery is expected to last 10 hours or 60 tests, whichever comes first, before needing to be recharged.
- The battery is expected to have a minimum life of 2 years under normal operating conditions (0° C to 45° C charging, -20° C to 60° C discharging).
- Store the battery in a cool dry place based on the following storage times at 65% ± 20% Relative humidity:
 - a. 1 month: -20° C to 60° C
 - b. 3 months: -20° C to 45° C
 - c. 12 months: -20° C to 25° C
- It is expected a full charge will take up to 12 hours.

Battery precautions

- **Note:** The battery cells can be shipped or stored down to -20° C but it is recommended to not be used until the temperature is at or greater than 0° C.
- It is recommended the battery is charged while not plugged in to the Cholestech LDX system.
- **Note:** Do not use the battery when the Cholestech LDX is connected to a printer.

 **Warning:**
HAZARDOUS LIVE PARTS INSIDE. DO NOT REMOVE COVER OR OPEN THE CASE. RISK OF ELECTRIC SHOCK. NO SERVICEABLE PARTS INSIDE.

 **Caution:**
DO NOT CONNECT SUPPLYING POWER AND LOAD DEVICE TO THE POWER BANK AT THE SAME TIME. DO NOT CONNECT THE OUTPUT PLUG OF THE POWER BANK TO ITS OWN INPUT JACK.

 **Important:**
DO NOT DISASSEMBLE BATTERY FOR ANY PURPOSE. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL. FOR BEST PERFORMANCE AND LONGEST PRODUCT LIFE, KEEP 40%~60% CHARGE STATE BEFORE STORAGE (LESS THAN 3 MONTHS). DO NOT USE THE PRODUCT IN TEMPERATURE BELOW 0° C OR ABOVE 60° C. AVOID DIRECT SUNLIGHT AND KEEP AWAY FROM FIRE. DO NOT IMMERSE THE PRODUCT IN WATER. KEEP OUT OF REACH OF CHILDREN.

Abbott
Cholestech LDX
Battery Kit

Quick Reference Guide

Size:
8.5 in. x 11 in.

Printed Colors



PMS 2925 C
Primary Blue



PMS 303 C
Dark Blue



Black



PMS 354 C
Green



PMS 185 C
Red

PN: 55170-4
Rev: AA

Date of Last Revision:
AA.4 2020/06/23