



Little Effort . . . Big Results

What leading dentists, hygienists, and scientists are saying about HydraBrush...

"The unique design eliminates many opportunities for user error and cleans the teeth more efficiently since all sides of both the upper and lower teeth are brushed simultaneously."

-Dr. Brad Fain, PhD
Director of the Department of Human Factors
Georgia Tech Research Institute

"Most people brush for 30-45 seconds and don't reach all areas. The HydraBrush makes what they do very effective no matter who is using the toothbrush, a child, an adult, someone with a disability. Whoever picks up that toothbrush will have the bristles in exactly the right spot every single time."

-Trisha O'Hehir, RDH
Editorial Director Hygienetown Magazine


"The HydraBrush is the most amazing preventative toothbrush advancement I have seen in my dental career."

-Howard Farran, DMD
Founder Dentaltown and Hygienetown

Patient compliance...
guaranteed!

ABRASION-PROOF

Short/quick (5mm) brush strokes
assure abrasion proof brushing as
confirmed by independent clinicals.



surrounds
teeth and gums

HYDRABRUSH IS THE BASS METHOD AUTOMATED

Certified ISO 13485 with Canadian requirements.

Big Results

- Automatically brushes according to the **Bass Brushing Technique** with 45° bristle angle along with short, quick brush strokes
- Teeth correctly position the brush, not the user, assuring patient compliance
- Brushes upper, lower, inside and outside surfaces simultaneously eliminating human error
- Gentle and total gum massage with every use. Dentists and Hygienists claim remarkable gingival improvement in as little as 3 weeks of HydraBrush use
- Longer, softer bristles contact gum tissue while shorter bristles clean the chewing surfaces
- 8 microbrushes completely surround all 6 tooth and gum surfaces at once, further assuring patient compliance
- Great for Orthodontic Patients. Effectively cleans around braces
- No timer needed. Brushing 40 seconds with HydraBrush is equivalent to several minutes with other power toothbrushes
- Clinically proven safe and effective. With 810 strokes per minute. HydraBrush provides effective interproximal cleaning. Visit our website for more information.
- Easier to work on patients (less blood/infection) after as little as 3 weeks of HydraBrush use
- Fresher breath, whiter teeth, better stain removal
- The only power toothbrush with the Arthritis Foundation "Ease-of-Use Commendation"
- Guaranteed patient compliance!



For Clinicals visit:
hydrabrush.com/clinicals

fast. effortless. precise.

Distributed by New Line Medical, Inc.

Call Now To Order

1-800-452-8909

Little Effort

- Just bite & guide
HydraBrush
for 40 seconds

Hygienetown

HydraBrush

The New Paradigm in Oral Hygiene




At last... patient compliance guaranteed!

As a dental professional, you know how difficult it can be to get your patients to brush properly, twice a day, 365 days a year. Time in the chair with them is not fun - red, swollen, inflamed gums cause pain and discomfort for the patient, mounting stress for you. Bleeding/tender gums make your job more difficult.

Help has arrived.

Patient compliance is now automatic and just *40 seconds away!*



hydra brush

Little Effort . . . Big Results

Hydrabrush Electric Toothbrush



HB001-108A-R
\$50.95
\$47.95 (6+)

Retail Price \$79.95

Brush Head Xtra Reach Standard Soft



HB005-EHEADER
\$7.00

Retail Price \$9.95

Brush Head Xtra Reach Ultra Soft



HB010-EHEADES
\$7.00

Retail Price \$9.95

HydraBrush 108A-R Model Includes

- 1 HydraBrush Power Handle
- 1 Countertop Pedestal with Hygienic Brush Storage
- 2 HydraBrush Brush heads (ultra soft + extra reach)
- 1 Compact charger and Cover
- 2 AA NiMH Rechargeable Batteries
- 1 Tongue Scraper
- 1 Travel Pouch

Brush Head Standard Soft



HB015-EHEADSS
\$7.00

Retail Price \$9.95

Brush Head Ultra Soft



HB020-EHEADUS
\$7.00

Retail Price \$9.95

*The ONLY toothbrush with
American Arthritis Foundation
Ease of Use Commendation.*

fast

delivers effective and convenient brushing in 40 seconds or less with 8 microbrushes surrounding teeth and massaging gums

effortless

unique technology simultaneously brushes upper, lower, inside and outside of tooth and gum surfaces - automatically!

precise

automatically brushes using the Bass Brushing Technique with angled bristles and short, quick brush strokes

HydraBrush, Inc.
701 S Andreason Dr., Suite C
Escondido, CA 92029



hydra brush

Little Effort . . . Big Results



smile

help has arrived

Certified ISO 13485 with Canadian requirements.

Distributed by New Line Medical, Inc. Call Now To Order 1-800-452-8909

"The HydraBrush" is the most amazing preventative dental toothbrush advancement I have ever seen in my dental career."

- Dr. Howard Farran, Founder & CEO
DentalTown.com & DentalTown Magazine

POPULAR MECHANICS

Electric toothbrushes have been buzzing, spinning and scouring our teeth in various ways for years, but the HydraBrush is the first to clean top and bottom teeth simultaneously. The rechargeable toothbrush uses two heads that surround and scrub your choppers while massaging your gums.

2005 FROST & SULLIVAN
Oral & Dental Care
Product of the Year --- HydraBrush