

Rapid, reliable results

Flexible with multiple  
sample options

Easy-to-use

Cost effective

## RELY™ H. pylori Rapid Test

Provides simple, fast and  
cost-effective detection  
of H. pylori, helping  
to lower the risk of  
developing a more  
serious illness



CLIA '88  
WAIVED



STANBIO  
LABORATORY

An EKF Diagnostics Company

## RELY™ H. pylori Rapid Test

# RELY™ H. pylori Rapid Test

After a brief and transient clinical illness, H. pylori infections are asymptomatic, yet clinical trials have demonstrated that H. pylori is associated with 90% of patients with duodenal ulcer and approximately 80% of those with gastric ulcer.<sup>1</sup>

### Enhanced sensitivity\*

87%

### Specificity\*

86%

### Reliable

Procedural and external controls provided

### Procedural steps

Two step test provides results in 10 minutes

### Sample

Whole blood (fingerstick or venipuncture), serum or plasma

### Easy-to-use

CLIA waived, requires minimal training

**Waived:** Whole blood

**Non-waived:** Serum and plasma

### Cost effective

Convenient room temperature/refrigerated storage and extended dating help manage inventory during seasonal shifts in testing

**Storage:** 2-30°C (36-86°F)

**Shelf life:** 18 Months from date of manufacture

### Easy and reliable detection of H. pylori means enhanced patient outcomes

- Minimizes the use of inappropriate treatment protocols
- Early detection and effective management of H. pylori reduce patient visits and cost of clinical follow-up
- Helps lower the risk of developing a more serious illness



Reference No.  
6300-020 (20 tests)



### Distributed by:

**New Line Medical, Inc.**

**800.452.8909**

**www.NewLineMedical.com**

**For over 50 years Stanbio Laboratory has been the leader in clinical diagnostics and is committed to provide only the highest quality products.**



CLIA '88  
WAIVED

CODE\*\*: 86318QW (Whole Blood)  
86318 (Serum/Plasma)



**STANBIO**  
LABORATORY

An EKF Diagnostics Company

ISO 13485:  
2003 Certified

RSS0015.03 Eff. Date: 6/8/12 DAR #: 12.0775 ©2012 Stanbio Laboratory


<sup>1</sup>CDC, Helicobacter pylori Fact Sheet for Health Care Providers. \*Whole blood, Fingerstick  
\*\*Exemplar code provided. Refer to regional CMS fee schedule for relevant codes and reimbursement rates.

# Helicobacter Pylori (H. pylori)

## What is it?

### Helicobacter Pylori

(H. Pylori) is a bacterium that causes ulcers in the stomach and duodenum as well as gastritis, or inflammation in the stomach.



**2/3**

world population infected

**#1**

Cause of ulcers

**25**

million Americans infected with H. pylori




**2x-6x**

chances of gastric cancer if infected

### Symptoms of H. pylori infection

- Belching
- Bloating
- Nausea, Vomiting
- Abdominal pain
- Fatigue
- Heartburn

Testing for H. Pylori is as easy as 1, 2, 3



- Collect whole blood, serum or plasma
- 2 drops of venipuncture whole blood. 1 drop of buffer
- 10 minutes to result

Distributed by:

New Line Medical, Inc.  
800.452.8909  
[www.NewLineMedical.com](http://www.NewLineMedical.com)

STANBIO<sup>®</sup>Chemistry

**EKF**  
Diagnostics for life