

**#1** Brand Recommended By Pharmacists and Doctors<sup>1</sup>

**OMRON**<sup>®</sup>


HEM-790IT

# Automatic Blood Pressure Monitor with ComFit™ Cuff



American Heart  
Association   
*Learn and Live™*



- **Detects Morning Hypertension**
  - Important Predictor of Increased Risk of Stroke
-  **IntelliSense™ Technology for Accurate, More Comfortable Readings**
- **Simple, Silent, 1-Touch Automatic Operation**
- **Up to 2.5 Times Faster Readings than Competitive Products<sup>2</sup>**
- **Extra Large Digital Display**
  - Systolic and Diastolic Blood Pressure with Pulse
- **Advanced Diagnostics**
  - Hypertension Indicator (According to the AHA guidelines)
  - Detects Morning Hypertension ▲
  - Detects Irregular Heartbeat ☹
  - 2 User - 200 Total Memory with Date and Time Stamp
    - 84 base memory per user
    - 8 weeks morning averages per user ☀
    - 8 weeks evening averages per user ☾
  - Two Modes
    - TruRead™ mode takes 3 consecutive measurements with AHA recommended 1 minute intervals and automatically displays the average. Intervals can be changes to 15, 30, 60 or 120 seconds.
    - Single Read Mode with Advanced Averaging takes individual readings and automatically displays the average of the last 3 readings taken within 10 minutes of each other.
- **ComFit™ Cuff with SigmaCuff Technology**
  - Pre-formed and Expandable to Easily Fit Regular and Large Arms
  - Fits Arms 9" to 17" in Circumference
  - Proprietary SigmaCuff Technology Allows Narrow Profile
    - Fast
    - Accurate
    - More Comfortable
- **AC Adapter, Carrying Case, Illustrated English/Spanish Instruction Manual, CD ROM, USB Cable and 4 "AA" Batteries Included**
- **5 Year Limited Warranty**

CLINICALLY PROVEN  
ACCURATE 

Automatic Blood Pressure Monitor



Om001-HEM-790IT

\$120.99/ bx.

New Line Medical, Inc.  
(800) 452-8909  
www.NewLineMedical.com